SlideShaft®

Push Button Adjustable Shaft Length
Neuro-stimulation Probe



- **All Neuro-stimulating Probe Shaft Sizes in One Probe (130 mm to 270 mm)**
- **Adjust Push Button Probe Shaft Length to Surgical Procedure**
- **Realtime Intraoperative Customization of Probe Shaft Length According to Patient Anatomy**

Ordering Information:

SlideShaft® Stimulator Probes

ITEM CODE	DESCRIPTION
NP001	SlideShaft®, Single Use Disposable Flush Tip Probe; 0.2 mm tip with push button adjustable shaft length 130 mm to 270 mm electrode; 1.9 m wire with DIN 42 802 connector
NP002	SlideShaft®, Single Use Disposable Ball Tip Probe; 2.3 mm ball with push button adjustable shaft length 130 mm to 270 mm electrode; 1.9 m wire with DIN 42 802 connector

NeuralPathways®, LLC
2674 East Main Street, Unit 432
Ventura, CA 93003 USA
Phone: **805.290.2710** • Fax: **877.330.1727**
www.neuralpathways.com
email: **sales@neuralpathways.com**

Manufactured by:
Technomed Europe
Amerikalaan 71
6199AE Maastricht Airport
The Netherlands

Manufactured Exclusively for:
Neural Pathways®, LLC
2674 East Main Street, Unit 432
Ventura, CA 93003 USA